

HOW TO DONATE TO ICARE

Donate Online From Anywhere



WeChat



AliPay



www.icarechina.org
(PayPal / Credit Cards)

Donating From Within China

For cash donations, please contact Zhang Yongliang (English Name: Aly) for details (+86-15352029585) (she speaks Chinese and English)

For direct wire transfer to our ICARE RMB bank account in China:

Our Bank Account No: 104000013624.

Account's Name: **NFIA (INTERNATIONAL
CONSULTANCY AND RESOURCE ENTERPRISE)**

Our Bank's Name: **BANK OF CHINA**

Bank Address: Gansu Province Lanzhou City Chengguan District Tianshui South Road No. 525 (KOHIA525)

Swift Code: BKCHCNBJ660

Office Address in English:

Chengguan District Hongyunrunyuan Building C6, Unit 2, Room 203
Lanzhou City, Gansu Province 730010

Our Office Address in Chinese: **KOHIABC6C2203 PostalCode:730010**

Donating From Outside of China

Companies and Individuals may donate by wire transfer from any country to our ICARE-Australia overseas fund as follows:

Name of account: ICARE Australia Overseas Aid Fund

Bank: Commonwealth Bank

Bank Address: Shop 67-68 224 Prospect Highway Seven Hills NSW 2147

SWIFT Code: CTBAAU2S

BSB Number: 062224

Account Number: 10119016

You may also donate directly to our Australia Office by sending a check to: ICARE Australia, First Floor, 15 Park Road Seven Hills NSW 2147.. Contact: Alex Genovese: 61-2-9838 7778.

Donating From Australia

You may donate directly to our Australia Office by sending a check to: ICARE Australia, First Floor, 15 Park Road, Seven Hills NSW 2147. Contact: Alex Genovese: 61-2-9838 7778.

Donating To Our USA 501(c)(3) Affiliate

(check, credit card or PayPal)

You may donate directly to BRIDGE TO THE NATIONS (BTTN), an affiliate of ICARE located in USA. BTTN is a tax-deductible 501(c)(3) organization (Tax Exemption No. 42-1623149). You may donate through the BTTN website at www.bttan.us (which allows donations internationally by credit card or PayPal). Or, you may send a check to BRIDGE TO THE NATIONS, P.O. Box 1829, Wake Forest, NC 27588 (USA).